



MOBILIZE PS 101 POLE STRETCHER

- The Mobilize Emergency Pole Stretchers with patient holding straps are lightweight, high-strength emergency stretchers suitable
 for first aid response, industrial sites and casualty rooms. This stretcher is characterized by its foldable design, which enhances portability and ease of storage when not in use.
- Widely utilized by emergency medical services (EMS), rescue teams, and healthcare professionals, the two-fold stretcher plays a critical role in providing rapid and secure patient transport.
- They are made with Stain & Liquid resistant PVC coated technical nylon fabric that are heat sealed making them resistant to blood, body fluids and mildew and easy to clean.
- All Mobilize Emergency Pole Stretchers feature lightweight aluminium poles with contoured, flexible hand grips for easy handling. These stretchers also hinges at the centre permitting the stretcher to dual fold for compact storage.
- The frame is engineered to support the weight of a patient while remaining collapsible for convenient storage and swift deployment. The stretcher surface is typically composed of a sturdy and easy-to-clean material, ensuring both patient comfort and hygiene.
- The "two-fold" aspect refers to the stretcher's ability to be folded at two distinct points, allowing it to becompacted into a smaller size for storage and transportation.
- This folding capability is crucial for emergency responders and medical professionals who often operate in confined spaces or need to quickly move and deploy equipment.



FEATURES:

- Durable, Stain & Liquid resistant PVC coated technical fabric resists stains and does not absorb bodily fluids.
- Folds lengthwise and crosswise for compact storage.
- Rugged spreader bars and legs that provide support and comfort.
- Nonsolid handgrips for firm control while handling.
- Heavy duty aluminum-alloy poles for easy handling.
- Medical Device (93/42/EC) Class 1.

PRODUCT SIZE

205 X 53 X 13 cm

PRODUCT SIZE - FOLDED

102 X 17 X 8

SELF WEIGHT

4.8 kg

LOAD LIMIT

159 kg

GROSS WEIGHT

5.5 kg

LIFT TYPE

Horizontal



Anti-Slip Handles



Body Strap(2 Nos.)



Thickened aluminium for heavy duty loading



Supportive feet



Fabric, Stain resistent and waterproof



J.THADHANI & CO.

New #12/ Old #28, Stringers Street, Chennai - 600001, Tamilnadu, India. **©** 044 - 4262 5223



www.thadhanisafety.com





MOBILIZE PS 101 POLE STRETCHER

USAGE:

- O To use the two-fold stretcher, emergency responders unfold the stretcher, lock it into position, and secure the patient onto its surface. Once the patient is safely situated, the stretcher can be lifted and carried by multiple responders using built-in handles.
- The stretcher's folding mechanism is user-friendly, enabling quick setup even in high-pressure situations.
- Two-fold stretchers are designed to accommodate a range of patient sizes and weights. The compact nature of the folded stretcher also makes it suitable for storage in ambulances, emergency vehicles or healthcare facilities with limited space.

Available in standard emergency safety orange with carry bag.

STEPS FOR USAGE:

STEP 1

Take out the stretcher from the carrying case.



STEP 2

Spread the stretcher.



STEP 3

Push the supportive bar to the end by hands or by foot.





